


HARTZINNKAPSELN MIT FLACHEM BODEN



Art.-Nr	Ø	Höhe	Vol.	Material	Verpackungs- einheit – Stk.
	mm	mm	ml		
15010	3,20	4,00	0,020	Sn/Ag97,5	250
15020	3,20	6,00	0,040	Sn/Ag97,5	250
15030	3,30	5,00	0,030	Sn/Ag97,5	250
15040	4,00	6,00	0,070	Sn/Ag97,5	250
15050	5,00	8,00	0,120	Sn/Ag97,5	250
15060	5,00	9,00	0,150	Sn/Ag97,5	250
15070	5,00	12,00	0,220	Sn/Ag97,5	250
15100	10,00	10,00		Sn/Ag97,5	100